

TC Forceps

03






Adson Delicate	02
Adson Micro	02,07
Adson Par Style	02
Adson	02,07,13
Adson-Brown	07
Baker Mod. (bonney)	10
Bonney	06
Cushing Micro	03
Cushing	03,08
Cushing-Potts	03
Delicate	09
Evans	05
Ferris-Smith	10
Gerald Micro	03
Gerald	03,08
Gillies	09
Leibinger	13
Leibinger-Debaquey	11,12
Leibinger-Williams	09
Mc Indoe	05
Obwegeser	05
Potts-Smith Micro	MICRO05
Potts-Smith	05
Rees Bayonet	06
Semkin Micro	08
Semkin	08
Skinner Kelicate	09
Skinner Micro	09
Skinner	09
Standard Delicate	04
Standard0	04,10
Thorlakson	13
Tickner	06
Uda Mod.	05
Wangensteen	06
Wheeler	13
Disectting Forceps	20,24
Tissue Forceps	25,30

WWW.MSCO.CO.UK



- Wide selection of patterns for use in all types of surgical procedures
- safe and secure grasping of finest tissue
- precisely machined and hand-crafted for longlasting performance
- gold-plated ends for easy identification


Serrations: 0.5 mm



REGULAR

This tooth size is ideal for use in standard pattern needle holders like Mayo-Hegar, Mathieu, Bozeman a.o. And in needle-pulling tissue forceps. For holding large needles with sutures sized 3-0 and larger.


Serrations: 0.4 mm



FINE

Best for use in needle holders such as Crile-Wood, DeBakey, Ryder, a.o. And in delicate needle pulling tissue forceps. For holding needles and sutures between 4-0 and 6-0.


Serrations: 0.65 mm



COARSE

This coarsely serrated insert is placed in the most powerful needle holder patterns and wire twisters and assures maximum holding ability when using heavy needles or surgical wire.

Serrations: 0.2 mm



MICRO

This micro size is extremely effective for use in microsurgical needle holders like Castroviejo, Mini-Ryder, Microvascular, Lemmon and other microoriented needle holders. For use with the most delicate needles and sutures from size 7-0 to 10-0. Excellent delicate needle holding ability without minimal damage to needle or suture.

SMOOTH

The smooth insert is designed to hold the smaller suture sizes without damage to the suture or needle. Its needle holding ability is somewhat limited, however. Used primarily in needle holders for plastic surgery and microsurgery.



Dissecting Forceps



Tissue Forceps




Atraumatic Forceps

WWW.MSSCO.CO.UK



Dissecting Forceps with tungsten carbide edges
Pinzas de diseccion con filos de carbure tungsteno
Pinces a dissection avec des lames en tungstene



Regular

BC-03-101
12 cm


Adson



Fine

BC-03-102
12 cm

Adson Par Style



Smooth

BC-03-103
12 cm

Adson Delicate



Regular

BC-03-104
15 cm


Adson



Fine

BC-03-105
15 cm

Adson Micro



Smooth

BC-03-106
15 cm

WWW.MSCO.CO.UK



Dissecting Forceps with tungsten carbide edges
Pinzas de disección con filos de carburo tungsteno
Pinces a dissection avec des lames en tungstène



1.8 mm



Regular

BC-03-107
18 cm

Cushing-Potts



1.8 mm



Regular

BC-03-108
18 cm

Cushing



Fine

BC-03-109
18 cm

Cushing Micro



1.8 mm



Fine

BC-03-110
18 cm
With lock

Cushing



1.8 mm



Regular

BC-03-111
18 cm

Cushing



Fine

BC-03-112
17.5 cm

Gerald



Fine

BC-03-113
17.5 cm

Gerald Micro



Dissecting Forceps with tungsten carbide edges
Pinzas de disección con filos de carburo tungsteno
Pincers a dissection avec des lames en tungstène





Dissecting Forceps with tungsten carbide edges
Pinzas de diseccion con filos de carbure tungsteno
Pincas a dissection avce des lames en tungstene



Regular

BC-03-122 15 cm

USA Mod.



Regular

BC-03-123 18 cm
BC-03-124 21 cm
BC-03-125 24.5 cm

Potts-smith



Fine

BC-03-126 21 cm

Potts-smith Micro



Regular

BC-03-127 25.5 cm

Evans



Fine

BC-03-128 15 cm

Mc Indoe



BC-03-129 16 cm
BC-03-130 18 cm

Obwegeser

WWW.MSSCO.CO.UK



Dissecting Forceps with tungsten carbide edges
Pinzas de disección con filos de carburo tungsteno
Pincas a dissection avec des lames en tungstène



Tickner



Regular

BC-03-131 14.5 cm
BC-03-132 19.5 cm



Bonney



Coarse

BC-03-133 17 cm



Regular


BC-03-134 15.5 cm
BC-03-135 23 cm
BC-03-136 25.5 cm

Wangensteen



Regular

BC-03-137 13.5 cm



Regular

BC-03-138 15.5 cm

Rees Bayonet



Regular

BC-03-139 19.5 cm
BC-03-140 23.5 cm

Rees Bayonet

WWW.MSSCO.CO.UK



Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pincers a dissection a griffes avec des lames en tungstene



Adson



Regular
BC-03-141 12 cm



Adson-Brown

3 x 4 Teeth



Regular
BC-03-142 12 cm



Adson Micro



BC-03-143 12 cm



Adson



Regular
BC-03-144 15 cm



Adson-Brown

3 x 4 Teeth



Regular
BC-03-145 15 cm



Fine
BC-03-146 15 cm



WWW.MSSCO.CO.UK



Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pincers a dissection a griffes avec des lames en tungstene



Semkin



Regular

BC-03-147 14.5 cm
BC-03-148 16 cm



Semkin Micro



Fine

BC-03-149 14.5 cm
BC-03-150 16 cm



Semkin Micro



Fine

BC-03-151 18 cm



Gerald



Regular

BC-03-152 17.5 cm



Gerald




Regular

BC-03-153 17.5 cm



Cushing



Regular

BC-03-154 18 cm

WWW.MSSCO.CO.UK



Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pincas a dissection a griffes avce des lames en tungstene





Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pincers a dissection a griffes avec des lames en tungstene



Coarse
BC-03-161 14.5 cm

Standard



Coarse
BC-03-162 16 cm

Standard



Coarse
BC-03-163 18 cm

Standard



Coarse
BC-03-164 17 cm

Baker Mod. (Bonney)



Coarse
BC-03-165 17 cm

Ferries-Smith





Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pincers a dissection a griffes avec des lames en tungstene

DEBAKEY-LEIBINGER Tissue Forceps with atraumatic serrations at the distal tip and finely cut Tungsten Carbide inserts placed directly behind.
The serration pattern cut in to the tips will grasp delicate tissues firmly yet with minimal trauma. The additional Tungsten Carbide insert allow the use of these forceps for controlled needle "pulling" as the needle is drawn through the tissue. Particularly use full for delicate vascular and cardiac surgery.

Atraumatic Serrations

1.5 mm



Fine

BC-03-166 15 cm
BC-03-167 19.5 cm
BC-03-168 24 cm

2.0 mm



Fine

BC-03-169 15 cm
BC-03-170 19.5 cm
BC-03-171 24 cm

2.7 mm



Fine

BC-03-172 15 cm
BC-03-173 19.5 cm
BC-03-174 24 cm

3.5 mm



Fine

BC-03-175 19.5 cm
BC-03-176 24 cm



Leibinger-DeBakey



Leibinger-DeBakey



Leibinger-DeBakey



Leibinger-DeBakey



Tissue Forceps with tungsten carbide edges
Pinzas quirurgicas con filos de carbure tungsteno
Pinces a dissection a griffes avec des lames en tungstene

Atraumatic Serrations
Handles are shapes
to allow comfortable
finger placement.

2.0 mm



Fine

Leibinger-DeBakey
BC-03-177 15 cm

2.0 mm



Fine

Leibinger-DeBakey
BC-03-178 19.5 cm

2.0 mm



Fine

Leibinger-DeBakey
BC-03-179 24 cm




Regular

BC-03-180 12.5 cm
BC-03-181 15 cm



Adson



Regular

BC-03-186 29.5 cm

Regular

BC-03-182 13.5 cm
BC-03-183 19.5 cm
BC-03-184 22.5 cm



Wheeler

Fine

BC-03-185 19.5 cm



Leibinger



Thorlakson



MEDICAL & SURGICAL SUPPLY UK LTD
WWW.MSSCO.CO.UK

501 Crompton Way
Bolton, Greater Manchester
United Kingdom
BL1 8UP

Telephone/WhatsApp
+447539 550563
+447549 292013

Email
info@mssco.co.uk
sales@mssco.co.uk